Every Cell Matters!

A busy Bio-Lab needs an efficient Cell Counter

Logos Biosystems' new LUNA-FX7™ Automated Cell Counter





Unmatched cell counting accuracy

- All-new optics
- Increased counting volume for the lowest CV per count
- Multichannel pipette-compatible 8-channel slides
- Fast and precise autofocus
- More robust and sophisticated counting algorithms
- Customizable cell-detection protocols

Optimized for bioprocess production applications

- Quality control and validation software
- Range of standard validation slides

21 CFR Part 11 ready

- User access and rights management
- Online data storage and control
- Encrypted electronic records

	NEW				
Compatible slides	LUNA™ 1-Channel Slides	LUNA™ 8-Channel Slides	LUNA™ 3-Channel Slides	LUNA™ Cell Counting Slides / PhotonSlides™	LUNA™ Reusable Slides
		THITTHE THE			ioni ioni ioni
Sample throughtput	1 sample	Up to 8 samples	Up to 3 samples	Up to 2 samples	1 sample
Sample loading volume	50 μL	10 μL / chamber	10 μL / chamber	10 μL / chamber	10 μL
Maximum analysis volume	5 μL	0.5 μL / chamber	1.3 μL / chamber	1.3 μL / chamber	1.3 µL

